

**Biogenic silicon by wet-alkaline digestion****Date: June 13 2001****Page: 1 of 4****References**

Krause G.L., C.L. Schelske and C.O. Davis, 1983. Comparison of three wet-alkaline methods of digestion of biogenic silicon in water. *Freshwater Biology* 13: 73-81.

For increase of sensitivity measure absorbances in cells longer than 1 cm – remember to modify concentrations of the standard and the control samples to be within the linear range for this longer cell.

**Principle**

Wet-alkaline digestion

The water sample is filtrated through a 0.45 µm membrane filter (Nucleopore polycarbonate, Millipore HA cellulose acetate or any other membrane filter).

The biogenic silicon in the suspended matter on the filter is extracted with sodium hydroxide at 85°C.

The extract is neutralised with sulphuric acid.

Determination of silicon in the sample extract according to Mwanza Method for Dissolved Reactive Silicon by the heteropolyblue method.

**Important**

All solutions are to be stored in polyethylene bottles and as far as possible also to be prepared in plastic equipment.

All solutions are to be prepared in water low in silicate.

**Reagents**

All reagents are to be stored in polyethylene bottles and as far as possible also to be prepared in plastic equipment.

All reagents are to be prepared in water low in silicate.

**Si-low water****Sodium hydroxide (NaOH) 0.5 M**

Dissolve 20 g of NaOH in Si-low distilled water and dilute to 1000 mL with distilled water.

Store in polyethylene bottle at room temperature.

**Biogenic silicon by wet-alkaline digestion****Date: June 13 2001****Page: 2 of 4****Sulphuric acid (H<sub>2</sub>SO<sub>4</sub>), 2.5 M**

Add with caution and while mixing 136 ml of concentrated H<sub>2</sub>SO<sub>4</sub> (density 1.84 g/mL) to about 700 ml of Si-low distilled water. Let the solution get room temperature and dilute to 1000 ml with water.

Store in polyethylene bottle at room temperature.

**Preservation and storage**

Perform analysis on fresh samples.

**Analysis**

All routines are as far as possible to be performed in plastic equipment.

Filter-Blank = Extraction of clean, dry filter as sample (2 replicates)

Reagent-Blank = Extraction without filter (2 replicates)

1. Filtrate a well-known volume of sample through a 0.45 µm membrane filter.
2. Place the filter in a plastic beaker.
3. Add 25 mL of 0.5 M NaOH pre-heated to 85°C.
4. Place in oven at 85°C for 15 minutes.
5. Cool to room temperature.
6. Add with caution and while mixing 2.50 ml 2.5 M H<sub>2</sub>SO<sub>4</sub>.
7. Transfer by means of Si-low distilled water the solution quantitatively to a 50 mL volumetric plastic flask.
8. Rinse beaker thoroughly with small portions of Si-low distilled water and transfer all rinse water to the volumetric flask as well.
9. Fill to 50 mL with Si-low distilled water.
10. If solution is unclear then either filtrate or let it settle.

Proceed with determination of dissolved reactive silicon by the heteropolyblue method according to Mwanza Method for Dissolved Reactive Silicon.

**Biogenic silicon by wet-alkaline digestion****Date: June 13 2001****Page: 3 of 4****Calculations**

Calculate the average absorbance<sub>Filter-Blank</sub>

Subtract the average absorbance<sub>Filter-Blank</sub> from filter samples.

Using blank corrected sample absorbance read the concentration in the sample extract from the calibration curve. Multiply by the dilution factor to obtain the sample extract concentration, A mg/L.

Calculate silicon content in 50 mL of sample extract, B mg Si:

$$B = \frac{\text{Sample extract concentration (A)} * 50}{1000} \text{ mg Si}$$

Since all silicon in the sample extract originates from the filtrated sample calculate the concentration of biogenic silicon in the sample, C mg/L PBSi:

$$C = \frac{B * 1000}{\text{ml sample filtrated}} \text{ mg/L Particulate Biogenic Si}$$

**Definitions****Si Dissolved reactive silicon**

Analysis of dissolved reactive silicon on original sample or on filtrated sample

**PBSi Particulate biogenic silicon**

Analysis of alkali-extractable silicon on particles from filtration of a well-known volume of sample (i.e. make analysis on the filter containing particles and extract filter-blanks containing the same type and amount of filter and reagent-blank without filter).

**Biogenic silicon by wet-alkaline digestion****Date: June 13 2001****Page: 4 of 4****Comments****Low concentration**

Filtrate a large sample volume to get a detectable amount of biogenic silicon on the filter.

Measure absorbances in long cells to increase sensitivity and modify standards for low concentrations.

**Contamination**

Since glass equipment contains silica all routines are as far as possible to be performed with plastic equipment.

All solutions are to be prepared in water low in silicate.

All solutions are to be stored in polyethylene bottles.